

NANO ICU-BED PREMIUM

operated electrically



ORDER INFORMATION

Item no.	Description	Size
94506	NANO ICU-bed PREMIUM	2280 x 965 mm

MATERIAL COMPOSITION

- carbon steel covered with antibacterial polyester-epoxy powder varnish
- chassis with ABS covers
- headend, footend and side panels made of plastic

DESCRIPTION

- Bed is supported on two tubular columns
- Bed segments are filled with radiolucent HPL
- With control system
- Wired control & monitoring
- Autorecovery function in back and thigh circulation segments
- Autocontour function
- CPR in backrest segment (one drive)
- 4 bumpers in the corners
- 4 hooks for liquids
- Battery backup
- Electrical and manual adjustment options (see below)
- DIN-rails on bed frame for additional accessories
- **Included accessories: IV pole, hand holder**

ADJUSTMENT OPTIONS

- **ELECTRICAL**
 - Backrest (0-70°) and cardiac chair position
 - Height adjustability (390 - 770mm)
 - (Reverse-)Trendelenburg position (0-15°)
- **MANUAL**
 - Footrest through raster

TECHNICAL DATA

- Central braking system, single castors 150 mm
- Permissible load: 250 kg
- 12 cm mattress with waterproof, vapour permeable cover
- Four-segment bed surface (three of them adjustable)
- Adjustability of the backrest: 0-70°
- Adjustability of the upper leg section: 0-40°
- Minimum height: 435 mm
- Maximum height: 805 mm

Stretching ca. (in%):	non-applicable
Sterilisation:	non-sterile
Washable:	Wipe disinfection*
CE-marking:	yes
Customs tariff number:	9402.9000
Durability:	non-applicable
Classification:	I
Disposable product:	no

LABELLING



*with commercially available detergents and disinfectants

Status: 04-2020