

# NEXO BED STANDARD

## Hospital bed, electric



### ORDER INFORMATION

| Item no. | Description       | Size          |
|----------|-------------------|---------------|
| 94601    | NEXO bed Standard | 2150 x 990 mm |

### MATERIAL COMPOSITION

- carbon steel covered with antibacterial polyester-epoxy powder varnish
- chassis with ABS covers
- side panels, head- and footend made of plastic

### DESCRIPTION

- with softdrop side rails (divided, full length side rails, made of plastic with inserts without controls)
- Bed segments are filled with radiolucent HPL
- With control system
- Wired control & monitoring
- Autorecovery function in back and thigh circulation segments
- Autocontour function
- CPR in backrest segment (one drive)
- 4 bumpers in the corners
- 8 hooks for liquids
- Electrical and manual adjustment options (see below)

### ADJUSTMENT OPTIONS

- ELECTRICAL
  - Backrest (0-70°) and cardiac chair position
  - Height adjustability
  - Trendelenburg position (0-16°)
  - Reverse-Trendelenburg position (0-17°)
- MANUAL
  - Footrest through raster

### TECHNICAL DATA

- Central braking system, single castors 150 mm
- Permissible load: 250 kg
- 12 cm mattress with waterproof, vapour permeable cover
- Four-segment bed surface (three of them adjustable)
- Adjustability of the backrest: 0-70°
- Adjustability of the upper leg section: 0-40°
- Minimum height: 410 mm
- Maximum height: 840 mm

|                        |                    |
|------------------------|--------------------|
| Stretching ca. (in%):  | non-applicable     |
| Sterilisation:         | non-sterile        |
| Washable:              | Wipe disinfection* |
| CE-marking:            | yes                |
| Customs tariff number: | 9402.9000          |
| Durability:            | non-applicable     |
| Classification:        | I                  |
| Disposable product:    | no                 |

### LABELLING



\*with commercially available detergents and disinfectants

Status: 04-2020